

## Anti-sarcoglycan Alpha SGCA Antibody

Catalog Number: A04756-1


### Overview

Product Name	Anti-sarcoglycan Alpha SGCA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-sarcoglycan Alpha SGCA Antibody catalog # A04756-1. Tested in WB,IP,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16586


### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IP: 1:50-1:100 IHC: 1:50-1:200

## Anti-sarcoglycan Alpha SGCA Antibody (A04756-1) Images



Western blot analysis of alpha sarcoglycan on rat heart cells lysates using anti-alpha sarcoglycan antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded mouse heart tissue using anti-alpha sarcoglycan antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-sarcoglycan Alpha SGCA Antibody

For Research Use Only. Not for use in diagnostic procedures.