

## Anti-NR2E1 Antibody

Catalog Number: A04767


### Overview

Product Name	Anti-NR2E1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NR2E1 Antibody catalog # A04767. Tested in WB,IF,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y466

### Technical Details

Immunogen	Recombinant full length Human NR2E1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IF: 1:50-1:100 IP: 1:50-1:100

## Anti-NR2E1 Antibody (A04767) Images



Western blot analysis of extracts of various cell lines, using NR2E1 antibody.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NR2E1 Antibody

For Research Use Only. Not for use in diagnostic procedures.