

Anti-NIFK (D229) Antibody

Catalog Number: A04797

Overview

Product Name	Anti-NIFK (D229) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NIFK (D229) Antibody catalog # A04797. Tested in IF, IHC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BYG3

Technical Details

Form	Liquid
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB: 1:500~1:1000 IHC: 1:50~1:200 IF: 1:50~1:200</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-NIFK (D229) Antibody (A04797) Images




Figure 2. Immunohistochemistry validation of NIFK using Anti-NIFK (D229) Antibody (A04797).

Immunohistochemistry (IHC) analyzes of NIFK (D229) pAb in paraffin-embedded human breast carcinoma tissue at 1:100. For more protocol information of IHC




Figure 1. Western blotting validation for Anti-NIFK (D229) Antibody A04797

Western blot (WB) analysis of NIFK (D229) polyclonal antibody at 1:500 dilution

Lane1:A549 whole cell lysate

Lane2:Raw264.7 whole cell lysate

Lane3:H9C2 whole cell lysate

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration




Figure 3. Western blotting validation for Anti-NIFK (D229) Antibody A04797

Western blot (WB) analysis of NIFK (D229) pAb at 1:500 dilution

Lane1:HCT116 whole cell lysate(40ug)

Lane2:SGC7901 whole cell lysate(40ug)

Lane3:CT26 whole cell lysate(40ug)

Lane4:PC12 whole cell lysate(40ug)

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NIFK (D229) Antibody