

## Anti-Phospho-NIFK (T234) Antibody

Catalog Number: A04797T234

### Overview


Product Name	Anti-Phospho-NIFK (T234) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Phospho-NIFK (T234) Antibody catalog # A04797T234. Tested in IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BYG3

### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human NIFK around the phosphorylation site of Thr234. AA range:200-249
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:40000

## Anti-Phospho-NIFK (T234) Antibody (A04797T234) Images

---



Immunohistochemistry validation of NIFK using Anti-Phospho-NIFK (T234) Antibody (A04797T234). Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°C)

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-NIFK (T234) Antibody

For Research Use Only. Not for use in diagnostic procedures.