

Anti-Tenascin-R/TNR Antibody Picoband®

Catalog Number: A04810-1

About TNR

Tenascin-R is a protein that in humans is encoded by the TNR gene. Tenascin-R(TNR) is an extracellular matrix protein expressed primarily in the central nervous system. It is a member of the tenascin(TN) gene family, which includes at least 3 genes in mammals: TNC(or hexabrachion), TNX(TNXB), and TNR. The genes are expressed in distinct tissues at different times during embryonic development and are present in adult tissues. TNR has been detected predominantly in the central nervous system and is localized around motor neurons and on motor axons in the spinal cord, cerebellum, hippocampus, and olfactory bulb. It is suggested that tenascin-R has a role in initiating the detachment of neuroblasts from tangential chains and in initiating radial migration of the cells.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-Tenascin-R/TNR Antibody Picoband®
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Tenascin-R/TNR Antibody Picoband® catalog # A04810-1. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q92752

Technical Details

Immunogen	E.coli-derived human Tenascin-R/TNR recombinant protein (Position: E60-R1007). Human TNR shares 94.2% and 93.6% amino acid (aa) sequence identity with mouse and rat TNR, respectively.
-----------	---

Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml, Mouse, Rat ELISA, 0.1-0.5 ug/ml, -

Anti-Tenascin-R/TNR Antibody Picoband® (A04810-1) Images



Western blot analysis of Tenascin-R/TNR using anti-Tenascin-R/TNR antibody (A04810-1). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: rat brain tissue lysates, Lane 2: mouse brain tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-Tenascin-R/TNR antigen affinity purified polyclonal antibody (Catalog # A04810-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for Tenascin-R/TNR at approximately 180 kDa. The expected band size for Tenascin-R/TNR is at 150 kDa.

Anti-Tenascin-R/TNR Antibody

For Research Use Only. Not for use in diagnostic procedures.