

Anti-PSD93 DLG2 Antibody

Catalog Number: A04826-2


Overview

Product Name	Anti-PSD93 DLG2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PSD93 DLG2 Antibody catalog # A04826-2. Tested in WB,Flow Cytometry,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15700


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 FC: 1:50-1:100 IP: 1:10-1:50

Anti-PSD93 DLG2 Antibody (A04826-2) Images



Western blot analysis of PSD93 on different lysates using anti-PSD93 antibody at 1/1,000 dilution. Positive control:
Lane 1: HeLa Lane 2: Mouse Thymus Lane 3: SH-SY5Y



Flow cytometric analysis of SH-SY5Y cells with PSD93 antibody at 1/100 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PSD93 DLG2 Antibody

For Research Use Only. Not for use in diagnostic procedures.