

## Anti-GluR2/3 GRM3 Antibody

Catalog Number: A04827

### Overview


Product Name	Anti-GluR2/3 GRM3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GluR2/3 GRM3 Antibody catalog # A04827. Tested in IHC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14832

### Technical Details


Immunogen	Synthesized peptide derived from human GluR2/3.
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Immunogen affinity purified
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB, 1:500-1:3000 IHC, 1:50-1:100

For protocols please visit <https://www.bosterbio.com/protocol-and-troubleshooting/>

## Anti-GluR2/3 GRM3 Antibody (A04827) Images



Western blot analysis of extracts from mouse brain, using GluR2/3 antibody A04827.



Immunohistochemical analysis of paraffin-embedded human brain tissue using GluR2/3 antibody A04827.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GluR2/3 GRM3 Antibody