

Anti-CLP36/PDLIM1 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A04832-3-Fluoro647

About PDLIM1

Enables transcription coactivator activity. Acts upstream of or within regulation of transcription by RNA polymerase II. Part of transcription regulator complex. Is expressed in several structures, including alimentary system; branchial arch; genitourinary system; heart; and sensory organ. Orthologous to human PDLIM1 (PDZ and LIM domain 1).

Overview

Product Name	Anti-CLP36/PDLIM1 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O70400

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of mouse CLP36/PDLIM1. Mouse CLP36/PDLIM1 shares 90% and 95% amino acid (aa) sequence identity with human and rat CLP36/PDLIM1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CLP36/PDLIM1 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.