

Anti-GPCR LGR8/RXFP2 Antibody Picoband® Biotin Conjugated

Catalog Number: A04848-1-Biotin

About RXFP2

Relaxin/insulin-like family peptide receptor 2, also known as RXFP2, is a human G-protein coupled receptor. It is mapped to 13q13.1. This gene encodes a member of the GPCR (G protein-coupled, 7-transmembrane receptor) family. Mutations in this gene are associated with cryptorchidism. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-GPCR LGR8/RXFP2 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8WXD0

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human GPCR LGR8.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GPCR LGR8/RXFP2 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.