

## Anti-AR-Alpha 1B ADRA1B Antibody

Catalog Number: A04856-1


### Overview

Product Name	Anti-AR-Alpha 1B ADRA1B Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-AR-Alpha 1B ADRA1B Antibody catalog # A04856-1. Tested in WB,IF applications. This antibody reacts with Human,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35368

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 475-520 of Human AR-alpha1B.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IF: 1:50-1:200

## Anti-AR-Alpha 1B ADRA1B Antibody (A04856-1) Images



Western blot (WB) analysis of AR-alpha1B polyclonal antibody at 1:500 dilution. Lane 1: A549 cell lysate, Lane 2: HEK293T cell lysate, Lane 3: PC12 cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AR-Alpha 1B ADRA1B Antibody

For Research Use Only. Not for use in diagnostic procedures.