

Anti-Protein Wnt-2b Wnt2b Antibody

Catalog Number: A04879


Overview

Product Name	Anti-Protein Wnt-2b Wnt2b Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Protein Wnt-2b Wnt2b Antibody catalog # A04879. Tested in WB,IHC,ICC,IF,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q93097


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 ICC: 1:100-1:500 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-Protein Wnt-2b Wnt2b Antibody (A04879) Images



Western blot analysis of Wnt2b on SiHa cell lysates using anti-Wnt2b at 1/500 dilution.



ICC staining Wnt2b in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

2 Publications Citing This Product

1. PubMed ID: 33359712, An R,Wang X,Yang L,Zhang J,Wang N,Xu F,Hou Y,Zhang H,Zhang L. Polystyrene microplastics cause granulosa cells apoptosis and fibrosis in ovary through oxidative stress in rats.Toxicology.2020 Dec 24;449:152665. doi:10.1016/j.tox.2020.152665.Epub ahead of p
2. PubMed ID: 32463570, Xu T,Pan L,Li L,Hu S,Zhou H,Yang C,Yang J,Li H,Liu Y,Meng X,Li J.MicroRNA-708 modulates Hepatic Stellate Cells activation and enhances extracellular matrix accumulation via direct targeting TMEM88. J Cell Mol Med. 2020 Jul;24(13):7127-7140.doi: 10.1111/jc

Visit bosterbio.com/anti-protein-wnt-2b-wnt2b-antibody-a04879-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein Wnt-2b Wnt2b Antibody

For Research Use Only. Not for use in diagnostic procedures.