

## Anti-RGS16 (phospho-Y168) Antibody

Catalog Number: A04881Y168


### Overview

Product Name	Anti-RGS16 (phospho-Y168) Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-RGS16 (phospho-Y168) Antibody catalog # A04881Y168. Tested in WB applications. This antibody reacts with Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O15492

### Technical Details

Immunogen	Synthetic phosphopeptide derived from human RGS16 around the phosphorylation site of Tyr168.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-RGS16 (phospho-Y168) Antibody (A04881Y168) Images



Western blot (WB) analysis of p-RGS16 (Y168) pAb at 1:500 dilution  
Lane1: The Eye tissue lysate of Mouse(40µg)  
Lane2: The Eye tissue lysate of Rat(40µg)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RGS16 (phospho-Y168) Antibody

For Research Use Only. Not for use in diagnostic procedures.