

## Anti-CNOT7 (N19) Antibody

Catalog Number: A04884


### Overview

Product Name	Anti-CNOT7 (N19) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CNOT7 (N19) Antibody catalog # A04884. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UIV1

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 10-54 of Human CNOT7
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE)
Suggested Dilutions	WB: 1:500-1:1000

## Anti-CNOT7 (N19) Antibody (A04884) Images



Western blot (WB) analysis of CNOT7 (N19) Antibody at 1:1000 dilution Lane1:Myla2059 whole cell lysate(40ug) Lane2:HEK293T whole cell lysate(40ug) Lane3:Hela whole cell lysate(40ug) Lane4:PC12 whole cell lysate(40ug) Lane5:3T3-L1 whole cell lysate(40ug)

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CNOT7 (N19) Antibody

For Research Use Only. Not for use in diagnostic procedures.