

Anti-CLN5 Antibody

Catalog Number: A04893-2


Overview

Product Name	Anti-CLN5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CLN5 Antibody catalog # A04893-2. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75503


Technical Details

Immunogen	Recombinant fusion protein of human CLN5(NP_006484.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

Anti-CLN5 Antibody (A04893-2) Images



Western blot (WB) analysis of CLN5 pAb at 1:500 dilution
Lane1:PMVEC whole cell lysate(40ug) Lane2:AML-12 whole cell lysate(40ug)
Lane3:L02 whole cell lysate(40ug)
Lane4:A549 whole cell lysate(40ug)



Immunohistochemistry of paraffin-embedded human appendix using CLN5 antibody at dilution of 1:150 .Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CLN5 Antibody

For Research Use Only. Not for use in diagnostic procedures.