

Anti-MSK1 (phospho-S360) RPS6KA5 Antibody

Catalog Number: A04922S360


Overview

Product Name	Anti-MSK1 (phospho-S360) RPS6KA5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MSK1 (phospho-S360) RPS6KA5 Antibody catalog # A04922S360. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75582

Technical Details

Immunogen	Synthetic phosphopeptide derived from human MSK1 around the phosphorylation site of Serine 360.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-MSK1 (phospho-S360) RPS6KA5 Antibody (A04922S360) Images



Western blot (WB) analysis of p-MSK1 (S360) polyclonal antibody at 1:500 dilution Lane1:MCF-7 whole cell lysate Lane2:MCF-7 treated with LPS(100ng/ml) for 5 minutes whole cell lysate Lane3:MCF-7 treated with LPS(100ng/ml) for 30 minutes whole cell lysate Lane4:AML-12 whole cell lysate(40ug) Lane5:H9C2 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MSK1 (phospho-S360) RPS6KA5 Antibody

For Research Use Only. Not for use in diagnostic procedures.