

Anti-MSK1 (phospho S376) RPS6KA5 Antibody

Catalog Number: A04922S376-1


Overview

Product Name	Anti-MSK1 (phospho S376) RPS6KA5 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-MSK1 (phospho S376) RPS6KA5 Antibody catalog # A04922S376-1. Tested in WB,ICC,IHC,IP applications. This antibody reacts with Human,Rat.
Application	IP, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75582


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser376 of human MSK1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-MSK1 (phospho S376) RPS6KA5 Antibody (A04922S376-1) Images



ICC staining MSK1(phospho S376) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining MSK1(phospho S376) in LO2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MSK1 (phospho S376) RPS6KA5 Antibody

For Research Use Only. Not for use in diagnostic procedures.