

Anti-Smad9 Antibody Picoband®

Catalog Number: A04932-1

About Smad9

Mothers against decapentaplegic homolog 9 also known as SMAD9, SMAD8, and MADH6 is a protein that in humans is encoded by the SMAD9 gene. This gene encodes a member of a family of proteins that act as downstream effectors of the bone morphogenetic protein (BMP) signaling pathway. The encoded protein is phosphorylated by BMP receptors, which stimulates its binding to SMAD4 and translocation into the nucleus, where it functions as a regulator of transcription. Activity of this protein is important for embryonic development. Mutation of this gene results in defects in pulmonary vasculature.

Overview


Product Name	Anti-Smad9 Antibody Picoband®
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-Smad9 Antibody Picoband® catalog # A07303-1. Tested in WB, ELISA applications. This antibody reacts with Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9JIW5

Technical Details

Immunogen	E.coli-derived mouse Smad9 recombinant protein (Position: Q153-H235). Mouse Smad9 shares 59.2% and 96.4% amino acid (aa) sequence identity with human and rat Smad9, respectively.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Lyophilized

Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml, Mouse, Rat ELISA, 0.1-0.5 ug/ml, -

Anti-Smad9 Antibody Picoband® (A04932-1) Images



Western blot analysis of Smad9 using anti-Smad9 antibody (A04932-1). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: rat brain tissue lysates, Lane 2: rat C6 whole cell lysates, Lane 3: mouse brain tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-Smad9 antigen affinity purified polyclonal antibody (Catalog # A04932-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for Smad9 at approximately 55 kDa. The expected band size for Smad9 is at 49 kDa.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Smad9 Antibody

For Research Use Only. Not for use in diagnostic procedures.