

Anti-RGS14 Antibody

Catalog Number: A04976


Overview

Product Name	Anti-RGS14 Antibody
Reactive Species	Human
Description	Boster Bio Anti-RGS14 Antibody catalog # A04976. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43566

Technical Details

Immunogen	Synthetic peptide, corresponding Human RGS14.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-RGS14 Antibody (A04976) Images



Immunohistochemistry (IHC) analyzes of RGS14 polyclonal antibody in paraffin-embedded human liver cancer tissue at 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RGS14 Antibody

For Research Use Only. Not for use in diagnostic procedures.