

## Anti-SGMS1 Antibody

Catalog Number: A04981


### Overview

Product Name	Anti-SGMS1 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-SGMS1 Antibody catalog # A04981. Tested in WB,ICC/IF applications. This antibody reacts with Human,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86VZ5


### Technical Details

Immunogen	Recombinant fusion protein of human SGMS1(NP_671512.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:200

## Anti-SGMS1 Antibody (A04981) Images



Western blot analysis of extracts of A-549 cells, using SGMS1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 930s.



Immunofluorescence analysis of HeLa cells using SGMS1 Rabbit pAb at dilution of 1:100. Blue: DAPI for nuclear staining.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-SGMS1 Antibody

For Research Use Only. Not for use in diagnostic procedures.