

Anti-AIDA Antibody Picoband® HRP Conjugated

Catalog Number: A04989-1-HRP

About AIDA

Predicted to enable phosphatidylinositol binding activity. Acts upstream of or within negative regulation of JUN kinase activity. Predicted to be located in cytoplasm. Predicted to be active in membrane.

Overview

| | |
|----------------------|---|
| Product Name | Anti-AIDA Antibody Picoband® HRP Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | WB, IHC, ELISA |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | Q96BJ3 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human AIDA recombinant protein (Position: E54-R284). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AIDA Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.