

Anti-Laminin Alpha-4 (A511) LAMA4 Antibody

Catalog Number: A04999-1


Overview

Product Name	Anti-Laminin Alpha-4 (A511) LAMA4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Laminin Alpha-4 (A511) LAMA4 Antibody catalog # A04999-1. Tested in IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16363

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 480-530 of Human Laminin alpha-4.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200 IF: 1:50-1:200

Anti-Laminin Alpha-4 (A511) LAMA4 Antibody (A04999-1) Images



Immunohistochemistry (IHC) analyzes of Laminin alpha-4 (A511) pAb in paraffin-embedded human colorectal carcinoma tissue at 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Laminin Alpha-4 (A511) LAMA4 Antibody

For Research Use Only. Not for use in diagnostic procedures.