

Anti-PLAA Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A05010-1-Fluoro488

About PLAA

Predicted to enable ubiquitin binding activity. Involved in cellular response to lipopolysaccharide; macroautophagy; and positive regulation of phospholipase A2 activity. Located in cytoplasm; extracellular exosome; and nucleus.

Overview

Product Name	Anti-PLAA Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y263

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human PLAA.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PLAA Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.