

## Anti-SOX7 Antibody Picoband®

Catalog Number: A05013-1

### About SOX7

SOX7 belongs to SOX gene family and SOX proteins are transcription factors with critical roles in the regulation of diverse developmental processes. This gene is mapped to 8p23.1. SOX7 gene contains at least 2 exons. In cotransfected 293 cells, Sox7 reduced Wnt (see WNT1)/beta-catenin (see CTNNB1)-stimulated transcription. SOX7 is a potent activator of FGF3 transcription. It not only plays a role in the transcriptional regulation of genes expressed in the hemogenic endothelium but also blocks further differentiation into blood precursors. And it may be required for the survival of both hematopoietic and endothelial precursors during specification.

### Overview


Product Name	Anti-SOX7 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-SOX7 Antibody Picoband® catalog # A05013-1. Tested in ELISA, WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BT81

### Technical Details

Immunogen	E.coli-derived human SOX7 recombinant protein (Position: K57-D343).
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.

Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5ug/ml, Human ELISA, 0.1-0.5ug/ml, -

## Anti-SOX7 Antibody Picoband® (A05013-1) Images



Western blot analysis of SOX7 using anti-SOX7 antibody (A05013-1). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human Colo320 whole cell lysates, Lane 2: human PC-3 whole cell lysates, Lane 3: human A431 whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-SOX7 antigen affinity purified polyclonal antibody (Catalog # A05013-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for SOX7 at approximately 43 kDa. The expected band size for SOX7 is at 43 kDa.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-SOX7 Antibody

For Research Use Only. Not for use in diagnostic procedures.