

## Anti-SOX7 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A05013-1-Fluoro594

### About SOX7

SOX7 belongs to SOX gene family and SOX proteins are transcription factors with critical roles in the regulation of diverse developmental processes. This gene is mapped to 8p23.1. SOX7 gene contains at least 2 exons. In cotransfected 293 cells, Sox7 reduced Wnt (see WNT1)/beta-catenin (see CTNNB1)-stimulated transcription. SOX7 is a potent activator of FGF3 transcription. It not only plays a role in the transcriptional regulation of genes expressed in the hemogenic endothelium but also blocks further differentiation into blood precursors. And it may be required for the survival of both hematopoietic and endothelial precursors during specification.

### Overview

Product Name	Anti-SOX7 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9BT81

### Technical Details

Immunogen	E.coli-derived human SOX7 recombinant protein (Position: K57-D343).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SOX7 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.