

Anti-Thrombospondin 4 THBS4 Antibody

Catalog Number: A05024


Overview

Product Name	Anti-Thrombospondin 4 THBS4 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Thrombospondin 4 THBS4 Antibody catalog # A05024. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35443


Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 551-600.
Isotype	IgG
Form	Liquid
Purification	Immunogen affinity purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used IHC-P, 1:50-300</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-Thrombospondin 4 THBS4 Antibody (A05024) Images



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Thrombospondin 4 THBS4 Antibody