

Anti-FGL1 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A05044-2-Fluoro594

About FGL1

Fibrinogen-like 1 is a member of the fibrinogen family. This protein is homologous to the carboxy terminus of the fibrinogen beta- and gamma-subunits which contains the four conserved cysteines of fibrinogens and fibrinogen related proteins. However, this protein lacks the platelet-binding site, cross-linking region and a thrombin-sensitive site which are necessary for fibrin clot formation. This protein may play a role in the development of hepatocellular carcinomas. Four alternatively spliced transcript variants encoding the same protein exist for this gene.

Overview

| | |
|----------------------|---|
| Product Name | Anti-FGL1 Antibody Picoband® Fluoro594 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q08830 |

Technical Details

| | |
|---------------------|---|
| Immunogen | A synthetic peptide corresponding to a sequence in the middle region of human FGL1. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGL1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.