

Anti-Adducin gamma ADD3 Antibody

Catalog Number: A05046

About ADD3

May play an important role in cell fate determination during lung development and in spermatogenesis.

Li C.-S., Chae S.-C., Lee J.-H., Zhang Q., Chung H.-T.J. Hum. Genet. 51:292-297(2006)

Overview

Product Name	Anti-Adducin gamma ADD3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Adducin gamma ADD3 Antibody catalog # A05046. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal 5G12G2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UEY8

Technical Details


Immunogen	Recombinant protein within human USP28 aa 550-700.
Predicted Reactive Species	Bovine, Canine, Equine, Pig, Rabbit
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL.
Purification	Immunogen affinity purified
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

	Boster Bio's internal QC testing used: WB, 1:500-1:2000
--	--

Anti-Adducin gamma ADD3 Antibody (A05046) Images



Western Blot (WB) analysis of HepG2 cells using Adducin gamma Polyclonal antibody.



Western Blot (WB) analysis of specific cells using Adducin gamma Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Adducin gamma ADD3 Antibody