

## Anti-Transcription factor E2F6 E2F6 Antibody

Catalog Number: A05062


### Overview

Product Name	Anti-Transcription factor E2F6 E2F6 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Transcription factor E2F6 E2F6 Antibody catalog # A05062. Tested in WB,IHC,IF,IP,CHIP applications. This antibody reacts with Human,Mouse.
Application	ChIP, IP, IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75461


### Technical Details

Immunogen	Recombinant full length Human E2F6.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC/IF: 1:50-1:200 IP/CHIP: 1:50-1:200

## Anti-Transcription factor E2F6 E2F6 Antibody (A05062) Images



Western blot analysis of extracts of various cell lines, using E2F6 antibody.



Immunofluorescence analysis of HeLa cell using E2F6 antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transcription factor E2F6 E2F6 Antibody

For Research Use Only. Not for use in diagnostic procedures.