

Anti-JIP-1 MAPK8IP1 Antibody

Catalog Number: A05068


Overview

Product Name	Anti-JIP-1 MAPK8IP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-JIP-1 MAPK8IP1 Antibody catalog # A05068. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UQF2


Technical Details

Immunogen	Synthesized peptide derived from human JIP-1 around the non-phosphorylation site of T103.
Isotype	IgG
Form	Liquid
Purification	Immunogen affinity purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB, 1:500-1:2000</p> <p>For protocols please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-JIP-1 MAPK8IP1 Antibody (A05068) Images



Western Blot (WB) analysis of COLO205 cells using JIP-1 Polyclonal antibody.



Western Blot (WB) analysis of specific cells using JIP-1 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-JIP-1 MAPK8IP1 Antibody