

Anti-OLAH Antibody Picoband® Biotin Conjugated

Catalog Number: A05069-3-Biotin

About OLAH

Enables dodecanoyl-[acyl-carrier-protein] hydrolase activity; myristoyl-[acyl-carrier-protein] hydrolase activity; and palmitoyl-[acyl-carrier-protein] hydrolase activity. Involved in medium-chain fatty acid biosynthetic process. Predicted to be located in cytosol.

Overview

Product Name	Anti-OLAH Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NV23

Technical Details

Immunogen	E.coli-derived human OLAH recombinant protein (Position: E13-F265).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-OLAH Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.