

## Anti-RAI3/GPRC5A Antibody Picoband® HRP Conjugated

Catalog Number: A05072-2-HRP

### About GPRC5A

This gene encodes a member of the type 3 G protein-coupling receptor family, characterized by the signature 7-transmembrane domain motif. The encoded protein may be involved in interaction between retinoid acid and G protein signalling pathways. Retinoic acid plays a critical role in development, cellular growth, and differentiation. This gene may play a role in embryonic development and epithelial cell differentiation.

### Overview

Product Name	Anti-RAI3/GPRC5A Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8NFJ5

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human RAI3/GPRC5A. Human RAI3/GPRC5A shares 95.2% amino acid (aa) sequence identity with mouse RAI3/GPRC5A.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAI3/GPRC5A Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.