

Anti-RTCB Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A05076-2-Fluoro488

About RTCB

RNA 2',3'-cyclic phosphate and 5'-OH ligase is a protein that in humans is encoded by the RTCB gene. Enables RNA ligase (ATP) activity and vinculin binding activity. Involved in tRNA splicing, via endonucleolytic cleavage and ligation. Located in cytosol and nucleoplasm. Part of tRNA-splicing ligase complex.

Overview

Product Name	Anti-RTCB Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y3I0

Technical Details

Immunogen	E.coli-derived human RTCB recombinant protein (Position: H227-G505).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RTCB Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.