

Anti-CD180/Rp105 Antibody

Catalog Number: A05085


Overview

| | |
|----------------------|--|
| Product Name | Anti-CD180/Rp105 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-CD180/Rp105 Antibody catalog # A05085. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat. |
| Application | ELISA, WB |
| Clonality | Polyclonal |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q99467 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Synthesized peptide derived from CD180 antigen at AA range: 221-270 |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Immunogen affinity purified |
| Suggested Dilutions | WB 1:500-1:2000 ELISA 1:10000 |

Anti-CD180/Rp105 Antibody (A05085) Images



Western Blot analysis of mouse brain cells using CD180 Polyclonal Antibody. Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:20000



Western Blot analysis of MOUSE-BRAIN cells using CD180 Polyclonal Antibody diluted at 1:1000. Secondary antibody was diluted at 1:20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD180/Rp105 Antibody

For Research Use Only. Not for use in diagnostic procedures.