

## Anti-SOAT1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A05096-1-Fluoro550

### About SOAT1

Sterol O-acyltransferase (acyl-Coenzyme A: cholesterol acyltransferase) 1, also known as SOAT1, is an enzyme that in humans is encoded by the SOAT1 gene. The protein encoded by this gene belongs to the acyltransferase family. It is located in the endoplasmic reticulum, and catalyzes the formation of fatty acid-cholesterol esters. This gene has been implicated in the formation of beta-amyloid and atherosclerotic plaques by controlling the equilibrium between free cholesterol and cytoplasmic cholesteryl esters. Alternatively spliced transcript variants have been found for this gene.

### Overview

Product Name	Anti-SOAT1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P35610

### Technical Details

Immunogen	E.coli-derived human SOAT1 recombinant protein (Position: R41-Y548). Human SOAT1 shares 88.2% and 87.4% amino acid (aa) sequence identity with mouse and rat SOAT1, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SOAT1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.