

Anti-DNA pol Delta 3 POLD3 Antibody

Catalog Number: A05111-1


Overview

Product Name	Anti-DNA pol Delta 3 POLD3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DNA pol Delta 3 POLD3 Antibody catalog # A05111-1. Tested in WB,ICC applications. This antibody reacts with Human,Mouse,Rat.
Application	ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15054

Technical Details

Immunogen	Synthetic peptide, corresponding to the 421-466 of Human DNA pol delta 3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 ICC: 1:50-1:200

Anti-DNA pol Delta 3 POLD3 Antibody (A05111-1) Images



DNA pol δ 3

Western blot (WB) analysis of DNA pol delta 3 pAb at 1:500 dilution Lane1:HuT78 whole cell lysate(40ug) Lane2:Myla2059 whole cell lysate(40ug) Lane3:SGC7901 whole cell lysate(40ug) Lane4:The Testis tissue lysate of Mouse(40ug) Lane5:PC12 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNA pol Delta 3 POLD3 Antibody

For Research Use Only. Not for use in diagnostic procedures.