

Anti-Citrate synthetase CS Antibody

Catalog Number: A05128-1


Overview

Product Name	Anti-Citrate synthetase CS Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Citrate synthetase CS Antibody catalog # A05128-1. Tested in WB,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75390


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2000 IF: 1:50-1:200

Anti-Citrate synthetase CS Antibody (A05128-1) Images



Western blot analysis of Citrate synthetase on Hela cell using anti-Citrate synthetase antibody at 1/1,000 dilution. Positive control: Lane 1: Rat kidney tissue Lane 2: Hela Lane 3: Mouse brain tissue



Immunohistochemical analysis of paraffin-embedded rat epididymis tissue using anti-Citrate synthetase antibody. Counter stained with hematoxylin.

1 Publications Citing This Product

1. PubMed ID: 32801864, Cai Z,Deng Y,Ye J,Zhuo Y,Liu Z,Liang Y,Zhang H,Zhu X,Luo Y,Feng Y,Liu R,Chen G,Wu Y,Han Z,Liang Y,Jiang F,Zhong W.Aberrant Expression of Citrate Synthase is Linked to Disease Progression and Clinical Outcome in Prostate Cancer.Cancer Manag Res.2020 Jul 22

Visit bosterbio.com/anti-citrate-synthetase-cs-antibody-a05128-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Citrate synthetase CS Antibody

For Research Use Only. Not for use in diagnostic procedures.