

## Anti-AR-Alpha 2C (R369) ADRA2C Antibody

Catalog Number: A05170-2

### About ADRA2C

Phosphoinositide-binding protein required for multivesicular body formation. Specifically binds phosphatidylinositol 3-phosphate (PtdIns(P3)). Plays a role in protein transport between cellular compartments. Promotes stability and cell surface expression of epithelial sodium channel (ENAC) subunits SCNN1A and SCNN1G By similarity. Not involved in EGFR degradation.

Haft C.R., Mol. Cell. Biol. 18:7278-7287(1998).

Xu Y., Nat. Cell Biol. 3:658-666(2001).

Hayama A., Submitted (AUG-2000) to the EMBL/GenBank/DDBJ databases.

### Overview

Product Name	Anti-AR-Alpha 2C (R369) ADRA2C Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-AR-Alpha 2C (R369) ADRA2C Antibody catalog # A05170-2. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P18825

### Technical Details

Immunogen	Synthesized peptide derived from internal of human SNX3.
Predicted Reactive Species	Boar, Bovine, Canine, Golden Hamster
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).

**Suggested Dilutions**

Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.

Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:

WB: 1:500-1:1000

IHC: 1:50-1:200

IF: 1:50-1:200

**Submit a product review to Biocompare.com**

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AR-Alpha 2C (R369) ADRA2C Antibody