

## Anti-IgD IGHD Antibody

Catalog Number: A05183

### Overview

Product Name	Anti-IgD IGHD Antibody
Reactive Species	Human
Description	Boster Bio Anti-IgD IGHD Antibody catalog # A05183. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01880

### Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 111-160.
Isotype	IgG
Form	Liquid
Purification	Immunogen affinity purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used IHC-P, 1:50-300</p> <p>For protocols please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a></p>

## Anti-IgD IGHD Antibody (A05183) Images

---



Immunohistochemical analysis of paraffin-embedded human tonsil, antibody was diluted at 1:200.

## Submit a product review to [Biocompare.com](http://Biocompare.com)

---

Submit a review of this product to [Biocompare.com](http://Biocompare.com) to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IgD IGHD Antibody