

Anti-SLC30A1 Antibody Picoband® HRP Conjugated

Catalog Number: A05201-1-HRP

About SLC30A1

Predicted to enable calcium channel inhibitor activity and zinc:proton antiporter activity. Involved in defense response to bacterium; intracellular zinc ion homeostasis; and zinc ion import into organelle. Located in several cellular components, including Golgi apparatus; basolateral plasma membrane; and nuclear membrane.

Overview

Product Name	Anti-SLC30A1 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y6M5

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human SLC30A1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC30A1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.