

## Anti-Kv1.4 Kcna4 Antibody

Catalog Number: A05208-1


### Overview

Product Name	Anti-Kv1.4 Kcna4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Kv1.4 Kcna4 Antibody catalog # A05208-1. Tested in WB,IHC,ICC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q61423


### Technical Details

Immunogen	Peptide
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Peptide affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:100 ICC: 1:50-1:200 FC: 1:50-1:100

## Anti-Kv1.4 Kcna4 Antibody (A05208-1) Images



Western blot analysis of Kv1.4 on rat brain (1) and mouse brain (2) tissue lysate using anti-Kv1.4 antibody at 1/1,000 dilution.



ICC staining Kv1.4 in N2A cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Kv1.4 Kcna4 Antibody

For Research Use Only. Not for use in diagnostic procedures.