

## Anti-FKBP10 Antibody Picoband® APC Conjugated

Catalog Number: A05217-1-APC

### About FKBP10

FK506-binding protein 10 is a protein that in humans is encoded by the FKBP10 gene. The protein encoded by this gene belongs to the FKBP-type peptidyl-prolyl cis/trans isomerase (PPIase) family. This protein localizes to the endoplasmic reticulum and acts as a molecular chaperone. Alternatively spliced variants encoding different isoforms have been reported, but their biological validity has not been determined.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-FKBP10 Antibody Picoband® APC Conjugated  |
| Reactive Species     | Human  |
| Application          | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .   |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.   |
| Host                 | Rabbit   |
| Uniprot ID           | Q96AY3   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | E.coli-derived human FKBP10 recombinant protein (Position: R29-E532).   |
| Cross Reactivity    | No cross-reactivity with other proteins.                                |
| Isotype             | Rabbit IgG  |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | APC<br>Excitation Wavelength: 633-647 nm<br>Emission Wavelength: 660 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users.                    |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FKBP10 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.