

Anti-ACAD9 Antibody

Catalog Number: A05223


Overview

Product Name	Anti-ACAD9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ACAD9 Antibody catalog # A05223. Tested in WB,ELISA applications. This antibody reacts with Human,Mouse,Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H845

Technical Details

Immunogen	A synthesized peptide derived from Human ACAD9(Accession Q9H845), corresponding to amino acid residues S358-F385.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-3000 ELISA(peptide): 1:20000-1:40000

Anti-ACAD9 Antibody (A05223) Images



Western blot analysis of extracts of various cell lines, using ACAD9 antibody at 1:1000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ACAD9 Antibody

For Research Use Only. Not for use in diagnostic procedures.