

## Anti-Cbx8 Antibody Picoband® FITC Conjugated

Catalog Number: A05234-2-FITC

### About CBX8

CBX8 functions as a transcriptional repressor and has a role in DNA damage response. This gene is mapped to chromosome 17q25.3 based on an alignment of the CBX8 sequence with the genomic sequence (GRCh38).

### Overview

Product Name	Anti-Cbx8 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9HC52

### Technical Details

Immunogen	E.coli-derived human Cbx8 recombinant protein (Position: M1-R389).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cbx8 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.