

Anti-Calsequestrin 1 CASQ1 Antibody

Catalog Number: A05235


Overview

Product Name	Anti-Calsequestrin 1 CASQ1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Calsequestrin 1 CASQ1 Antibody catalog # A05235. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P31415


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-Calsequestrin 1 CASQ1 Antibody (A05235) Images



Immunohistochemical analysis of paraffin-embedded rat skeletal muscle tissue using anti-Calsequestrin 1 antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human fetal skeletal muscle tissue using anti-Calsequestrin 1 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calsequestrin 1 CASQ1 Antibody

For Research Use Only. Not for use in diagnostic procedures.