

## Anti-IL-31R alpha Antibody

Catalog Number: A05256


### Overview

Product Name	Anti-IL-31R alpha Antibody
Reactive Species	Human
Description	Boster Bio Anti-IL-31R alpha Antibody catalog # A05256. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal 4A11A6G11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NI17

### Technical Details

Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human COQ2
Predicted Reactive Species	Canine, Equine, Goat, Pig, Rabbit, Yeast, Zebrafish
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB, 1:500-1:2000</p>

## Anti-IL-31R alpha Antibody (A05256) Images



Western Blot (WB) analysis of extracts from Jurkat cells, using IL-31R Polyclonal antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-31R alpha Antibody