

## Anti-Calmegin/Clgn Antibody Picoband® Biotin Conjugated

Catalog Number: A05261-2-Biotin

### About Clgn

Calmegin, also known as CLGN, is a protein which in humans is encoded by the CLGN gene. This gene belongs to the calreticulin family, which includes calreticulin, calnexin, and calmegin, and encodes a calcium-binding molecular chaperone specifically expressed in pachytene stage male germ cells. It is required for the proper folding of newly synthesized membrane proteins in the endoplasmic reticulum including those critical for sperm migration from the uterus into the oviduct and sperm adhesion to and penetration of the zona pellucida. This gene plays a key role in spermatogenesis and male infertility. Alternative splice variants exist for this gene.

### Overview

Product Name	Anti-Calmegin/Clgn Antibody Picoband® Biotin Conjugated
Reactive Species	Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P52194

### Technical Details

Immunogen	E.coli-derived mouse Calmegin/Clgn recombinant protein (Position: N18-D611).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calmegin/Clgn Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.