

Anti-Calmegin/Clgn Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A05261-2-Fluoro594

About Clgn

Calmegin, also known as CLGN, is a protein which in humans is encoded by the CLGN gene. This gene belongs to the calreticulin family, which includes calreticulin, calnexin, and calmegin, and encodes a calcium-binding molecular chaperone specifically expressed in pachytene stage male germ cells. It is required for the proper folding of newly synthesized membrane proteins in the endoplasmic reticulum including those critical for sperm migration from the uterus into the oviduct and sperm adhesion to and penetration of the zona pellucida. This gene plays a key role in spermatogenesis and male infertility. Alternative splice variants exist for this gene.

Overview

Product Name	Anti-Calmegin/Clgn Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P52194

Technical Details

Immunogen	E.coli-derived mouse Calmegin/Clgn recombinant protein (Position: N18-D611).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calmegin/Clgn Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.