

Anti-FLVCR1 Antibody

Catalog Number: A05270-3

About FLVCR1

FLVCR1(feline leukemia virus subgroup C receptor 1) has been reported to have a crucial role in a variety of biological processes, including cell proliferation, cell death, apoptosis, oxidative stress response, cellular differentiation, and metabolism (PMID: 29532854). In addition, FLVCR1 can be applied as a clinical prognostic marker for patients with ESCC (PMID: 33842377).

Overview

| | |
|----------------------|--|
| Product Name | Anti-FLVCR1 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-FLVCR1 Antibody catalog # A05270-3. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat. |
| Application | ELISA, IP, IF, IHC, ICC, WB |
| Clonality | Polyclonal |
| Formulation | 500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | Q9Y5Y0 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human FLVCR1 recombinant protein (Position: R3-P105). |
| Form | Liquid |
| Concentration | 500 ug/ml |
| Purification | Immunogen affinity purified. |
| Suggested Dilutions | Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000 |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FLVCR1 Antibody

For Research Use Only. Not for use in diagnostic procedures.