

Anti-Sarcoglycan delta SGCD Antibody

Catalog Number: A05297-1


Overview

Product Name	Anti-Sarcoglycan delta SGCD Antibody
Reactive Species	Human
Description	Boster Bio Anti-Sarcoglycan delta SGCD Antibody catalog # A05297-1. Tested in WB,IP,IHC applications. This antibody reacts with Human.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92629


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 IP: 1:10-1:50

Anti-Sarcoglycan delta SGCD Antibody (A05297-1) Images



Western blot analysis of delta Sarcoglycan on human lung cells lysates using anti- delta Sarcoglycan antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue using anti- delta Sarcoglycan antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Sarcoglycan delta SGCD Antibody

For Research Use Only. Not for use in diagnostic procedures.