

## Anti-HMGCS1 Antibody

Catalog Number: A05313


### Overview

Product Name	Anti-HMGCS1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HMGCS1 Antibody catalog # A05313. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q01581


### Technical Details

Immunogen	Recombinant fusion protein of human HMGCS1(NP_002121.4).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:2000 ICC/IF: 1:50-1:200

## Anti-HMGCS1 Antibody (A05313) Images



Western blot analysis of HMGCS1 expression in HeLa (A), Jurkat (B), mouse liver (C), rat spinal cord (D) whole cell lysates.



Immunofluorescent analysis of HMGCS1 staining in HeLa cells.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HMGCS1 Antibody

For Research Use Only. Not for use in diagnostic procedures.